

Specification Sheet



Product Name: TCS 2210

Item No.: MED18931

Purity Specification: $\geq 98\%$

Molecular Formula :

CAS Number: 1201916-31-5 *Formula Weight:*


Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816