

Specification Sheet



Product Name: Toceranib

Item No.: MED19347

Purity Specification: $\geq 98\%$

Molecular Formula : C₂₂H₂₅N₄O₂

CAS Number: 356068-94-5 **Formula Weight:** 396.458


Storage:

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816