

Specification Sheet



Product Name: TP-0903

Item No.: MED19420

Purity Specification: 98%


Molecular Formula : C₂₄H₃₀CLN₇O₂S

CASNumber: 1341200-45-0 **FormulaWeight:** 516.059

Storage:

| Test | Specification | |
|--------|---------------|------|
| | Min. | Max. |
| PURITY | 98% | 100% |

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816