

Specification Sheet



Product Name: Tandutinib (MLN518) HCl(free base)

Item No.: MED19459

Purity Specification: $\geq 98\%$

Molecular Formula : C₃₁H₄₂N₆O₄

CASNumber: 387867-13-2 **FormulaWeight:** 562.703


Storage:

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816