

Specification Sheet



Product Name: PSI

Item No.: MED19540

Purity Specification: 98%

Molecular Formula : C₃₂H₅₀N₄O₈

CAS Number: 158442-41-2

Formula Weight: 618.761


Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816