

Specification Sheet



Product Name: PIK-III

Item No.: MED19591

Purity Specification: $\geq 98\%$

Molecular Formula : C₁₇H₁₇N₇

CAS Number: 1383716-40-2 **Formula Weight:** 319.364

Storage: -20 ° C

Test	Specification	
	Min.	Max.
PURITY	98%	100%

Structure :

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816