

Specification Sheet



Product Name: Fmoc-Dap(Dnp)-OH

Item No.: MED19731

Purity Specification: $\geq 98\%$

Molecular Formula : C₂₄H₂₀N₄O₈

CAS Number: 140430-54-2 *Formula Weight:* 492.438


Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816