

Specification Sheet



Product Name: Arachidonic Acid Alkyne(solution)

Item No.: MED21306

Purity Specification: $\geq 98\%$

Molecular Formula :

CASNumber: 1219038-32-0 **FormulaWeight:**

Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816