

Specification Sheet



Product Name: Hyperforin

Item No.: MED21376

Purity Specification: 98%

Molecular Formula : C₃₅H₅₂O₄

CAS Number: 11079-53-1

Formula Weight: 536.785


Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816