

Specification Sheet



Product Name: N-heptanoyl-L-Homoserine lactone

Item No.: MED21443

Purity Specification: 98%

Molecular Formula : C₁₁H₁₉NO₃


CASNumber: 177158-20-2

FormulaWeight: 213.273

Storage: -20 ° C

Test	Specification	
	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816