

Specification Sheet



Product Name: tetranor-PGAM(solution)

Item No.: MED21515

Purity Specification: 98%

Molecular Formula :

CAS Number: 52510-53-9

Formula Weight:

Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816