

Specification Sheet



Product Name: Digitoxigenin

Item No.: MED21787

Purity Specification: 98%

Molecular Formula : C₂₃H₃₄O₄


CAS Number: 143-62-4

Formula Weight: 374.514

Storage: -20 ° C

| Test | Specification | |
|--------|---------------|------|
| | Min. | Max. |
| PURITY | 98% | 100% |

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816