

# Specification Sheet



**Product Name:** Docosaheptaenoyl PAF C-16(solution)

*Item No.:* MED21933

**Purity Specification:** 98%

**Molecular Formula :** C<sub>36</sub>H<sub>69</sub>NO<sub>7</sub>P

**CASNumber:** 132213-85-5

**FormulaWeight:** 658.909


**Storage:** -20 ° C

## Test

## Specification

|        | Min. | Max. |
|--------|------|------|
| PURITY | 98%  | 100% |

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816