

Specification Sheet



Product Name: CRANAD 28

Item No.: MED21941

Purity Specification: $\geq 98\%$

Molecular Formula :

CASNumber: 1623747-97-6 *FormulaWeight:*

Storage: -20 ° C

Test	Specification	
	Min.	Max.
PURITY	98%	100%

Structure :

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816