

Specification Sheet



Product Name: ($\pm$)5-HETE lactone(solution)

Item No.: MED21998

Purity Specification: $\geq 98\%$

Molecular Formula :

CASNumber: 73279-37-5

FormulaWeight:

Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816