

Specification Sheet



Product Name: PTP Inhibitor III

Item No.: MED22067

Purity Specification: $\geq 95\%$

Molecular Formula :

CASNumber: 29936-81-0

FormulaWeight:

Storage: -20 ° C

Test

Specification

	Specification	
	Min.	Max.
PURITY	95%	100%

Structure :

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816