

Specification Sheet



Product Name: DA-3003-2

Item No.: MED22186

Purity Specification: $\geq 98\%$

Molecular Formula :

CAS Number: 383907-47-9 **Formula Weight:**

Storage: -20 ° C

Test	Specification	
	Min.	Max.
PURITY	98%	100%

Structure :

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816