

Specification Sheet



Product Name: 3,3'-Diiodo-L-thyronine

Item No.: MED22214

Purity Specification: $\geq 98\%$

Molecular Formula : C₁₅H₁₃I₂NO₄

CAS Number: 4604-41-5

Formula Weight: 525.077


Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816