

Specification Sheet



Product Name: 2-Picolinic acid

Item No.: MED22300

Purity Specification: 98%

Molecular Formula : C₆H₅NO₂

CAS Number: 98-98-6

Formula Weight: 123.109


Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816