

Specification Sheet



Product Name: D(-)-Arginine

Item No.: MED22384

Purity Specification: 98%

Molecular Formula : C₆H₁₄N₄O₂

CAS Number: 157-06-2

Formula Weight: 174.201


Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816