

Specification Sheet



Product Name: Isoproterenol sulfate dihydrate

Item No.: MED22568

Purity Specification: $\geq 98\%$

Molecular Formula : C₂₂H₄₀N₂O₁₂S

CASNumber: 299-95-6

FormulaWeight: 556.624


Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816