

Specification Sheet



Product Name: L-Arginine (hydrochloride)

Item No.: MED22600

Purity Specification: 98%

Molecular Formula : C₆H₁₅CLN₄O₂

CAS Number: 1119-34-2

Formula Weight: 210.662


Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816