

Specification Sheet



Product Name: Phenylalanine & tyrosine metabolism

Item No.: MED22704

Purity Specification: $\geq 98\%$

Molecular Formula : C₉H₈O₃

CASNumber: 156-06-9

FormulaWeight: 164.158


Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816