

Specification Sheet



Product Name: Vitamin D3-d7

Item No.: MED25358

Purity Specification: $\geq 90\%$

Molecular Formula : $C_{61}H_{99}N_2 \cdot ClO_4$

CASNumber: 1627523-19-6 **FormulaWeight:** 959.900

Storage: -80°C

Test

Specification

	Min.	Max.
PURITY	90%	100%

Structure :

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816