

Specification Sheet



Product Name: BAY-6035

Item No.: MED25603

Purity Specification: $\geq 95\%$

Molecular Formula : C₂₂H₂₈N₄O₃

CASNumber: 2247890-13-5 FormulaWeight: 396.480

Storage: -20°C

Test

Specification

	Min.	Max.
PURITY	95%	100%

Structure :

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816