

Specification Sheet



Product Name: Docosahexaenoic Acid-d5 Standard

Item No.: MED26221

Purity Specification: $\geq 98\%$

Molecular Formula : C₂₂H₃₂O₂

CASNumber: 1197205-71-2 **FormulaWeight:** 333.519


Storage: -20°C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816