

Specification Sheet



Product Name: 3-hydroxy-DL-Kynurenine

Item No.: MED26380

Purity Specification: 98%

Molecular Formula : C₁₀H₁₂N₂O₄

CAS Number: 484-78-6

Formula Weight: 224.213


Storage: -20

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816