

Specification Sheet



Product Name: 3-methylquinoline-8-sulfonic acid

Item No.: MED32390

Purity Specification: $\geq 99\%$

Molecular Formula :

CASNumber: 153886-69-2

FormulaWeight:

Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	99%	100%

Structure :

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816