

Specification Sheet



Product Name: 2-(p-tolyl)propanoic acid

Item No.: MED32811

Purity Specification: 100%

Molecular Formula : C₁₀H₁₂O₂

CAS Number: 938-94-3

Formula Weight: 164.201


Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	100%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816