

Specification Sheet



Product Name: 3-(m-tolyl)propanoic acid

Item No.: MED32960

Purity Specification: 99%

Molecular Formula : C₁₀H₁₂O₂

CAS Number: 3751-48-2

Formula Weight: 164.201


Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	99%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816