

# Specification Sheet



**Product Name:** pyridine-3-sulfonic acid

*Item No.:* MED33283

Purity Specification: 99%

Molecular Formula : C<sub>5</sub>H<sub>5</sub>NO<sub>3</sub>S


CASNumber: 636-73-7

FormulaWeight: 159.163

Storage: -20 ° C

Test	Specification	
	Min.	Max.
PURITY	99%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816