

# Specification Sheet



## Product Name:

Item No.: MED34143

Purity Specification: 100%

Molecular Formula : C<sub>18</sub>H<sub>20</sub>CLNOS

CASNumber: 910138-96-4

FormulaWeight: 333.876


Storage: -20 ° C

## Test

## Specification

|        | Min. | Max. |
|--------|------|------|
| PURITY | 100% | 100% |

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: [www.med-bio.cn](http://www.med-bio.cn) Tel: 400-186-1816