

Specification Sheet



Product Name: 14,15-dehydro Leukotriene B4

Item No.: MED50035

Purity Specification: 97%

Molecular Formula :

CASNumber: 114616-11-4

FormulaWeight:

Storage: -80 ° C

Test

Specification

	Min.	Max.
PURITY	97%	100%

Structure :

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816