

Specification Sheet



Product Name: AMG510

Item No.: MED50064

Purity Specification: $\geq 99\%$

Molecular Formula : C₃₀H₃₀F₂N₆O₃

CASNumber: 2296729-00-3 **FormulaWeight:** 560.590

Storage: -20 ° C

Test	Specification	
	Min.	Max.
PURITY	99%	100%

Structure :

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816