

Specification Sheet



Product Name: BIS-TRIS

Item No.: MED80036

Purity Specification: 99%

Molecular Formula : C₈H₁₉NO₅

CAS Number: 6976-37-0

Formula Weight: 209.240


Storage: Store at RT

Test

Specification

| | Min. | Max. |
|--------|------|------|
| PURITY | 99% | 100% |

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816