

Specification Sheet



Product Name: naringenin

Item No.: MED80064

Purity Specification: $\geq 98\%$

Molecular Formula : C₁₅H₁₂O₅

CASNumber: 480-41-1

FormulaWeight: 272.253


Storage: 2~8 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816