

Specification Sheet



Product Name: Icariin

Item No.: MED80093

Purity Specification: 98%

Molecular Formula : C₃₃H₄₀O₁₅

CAS Number: 489-32-7

Formula Weight: 676.662


Storage: 2~8 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816