

# Specification Sheet



**Product Name:** Polydatin

*Item No.:* MED80104

**Purity Specification:** 98%

**Molecular Formula :** C<sub>20</sub>H<sub>22</sub>O<sub>8</sub>

**CAS Number:** 65914-17-2

**Formula Weight:** 390.384


**Storage:** 2~8 °C

## Test

## Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: [www.med-bio.cn](http://www.med-bio.cn) Tel: 400-186-1816