

Specification Sheet



Product Name: Tetrandrine

Item No.: MED80114

Purity Specification: 98%

Molecular Formula : C₃₈H₄₂N₂O₆

CASNumber: 518-34-3

FormulaWeight: 622.750


Storage: 2~8 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816