

^1H NMR (400 MHz, CHCl_3 - d) δ ppm 0.93 (t, $J=7.40$ Hz, 3 H) 1.60 - 1.75 (m, 2 H) 2.23 (dtd, $J=12.45, 11.28, 11.28, 8.91$ Hz, 1 H) 2.52 (t, $J=7.28$ Hz, 2 H) 2.76 (dddd, $J=12.55, 8.72, 6.09, 1.51$ Hz, 1 H) 3.47 (s, 2 H) 4.28 (ddd, $J=10.98, 9.35, 6.02$ Hz, 1 H) 4.37 - 4.53 (m, 1 H) 4.59 (ddd, $J=11.42, 8.78, 6.65$ Hz, 1 H) 7.64 (br. s., 1 H)

