

^1H NMR (400 MHz, $\text{DMSO}-d_6$) δ ppm 12.11 (br. s., 1 H) 8.80 (s, 1 H) 8.68 (s, 1 H) 8.37 (s, 1 H) 7.50 - 7.78 (m, 1 H) 6.99 (dd, $J=3.51, 1.51$ Hz, 1 H) 4.54 (td, $J=9.66, 4.02$ Hz, 1 H) 3.24 - 3.32 (m, 1 H) 3.11 - 3.24 (m, 1 H) 2.36 - 2.49 (m, 1 H) 1.73 - 1.94 (m, 1 H) 1.49 - 1.73 (m, 3 H) 1.39 - 1.49 (m, 1 H) 1.25 - 1.39 (m, 2 H) 1.12 - 1.25 (m, 1 H)

