

MED50168-M2F0168AC-H-20220610




7.715
7.268
4.635
4.618
4.613
4.606
4.597
4.590
4.585
4.568
4.510
4.489
4.466
4.317
4.302
4.294
4.290
4.279
4.274
4.266
4.251
3.476
2.781
2.778
2.775
2.764
2.750
2.747
2.540
2.522
2.504
2.269
2.246
2.238
2.216
1.676
1.658
1.639
1.621
1.603
0.952
0.933
0.915
-0.000

Current Data Parameters
NAME MED50168-M2F0168AC
EXPNO 1
PROCNO 1

F2 - Acquisition Parameters
Date_ 20220610
Time 9.25
INSTRUM spect
PROBHD 5 mm PABBO BB/
PULPROG zg30
TD 65536
SOLVENT CDC13
NS 16
DS 2
SWH 8012.820 Hz
FIDRES 0.122266 Hz
AQ 4.0894465 sec
RG 102.33
DW 62.400 usec
DE 6.50 usec
TE 291.4 K
D1 1.00000000 sec
TD0 1

==== CHANNEL f1 =====
SFO1 400.1324710 MHz
NUC1 1H
P1 10.00 usec
PLW1 16.00000000 W

F2 - Processing parameters
SI 65536
SF 400.1300067 MHz
WDW EM
SSB 0
LB 0.30 Hz
GB 0
PC 1.00



0.91
1.03
0.96
1.02
1.84
1.03
1.83
1.03
3.13
3.00