


KDS-3-H-20211223



Current Data Parameters
NAME KDS-3-H-20211223
EXPNO 1
PROCNO 1

F2 - Acquisition Parameters
Date_ 20211223
Time 9.45
INSTRUM spect
PROBHD 5 mm PABBO BB/
PULPROG zg30
TD 65536
SOLVENT CDC13
NS 16
DS 2
SWH 8012.820 Hz
FIDRES 0.122266 Hz
AQ 4.0894465 sec
RG 50.71
DW 62.400 usec
DE 6.50 usec
TE 296.1 K
D1 1.0000000 sec
TD0 1

==== CHANNEL f1 =====
SFO1 400.1324710 MHz
NUC1 1H
P1 10.00 usec
PLW1 16.0000000 W

F2 - Processing parameters
SI 65536
SF 400.1300420 MHz
WDW EM
SSB 0
LB 0.30 Hz
GB 0
PC 1.00