

Product Name: Delipidized Serum, (Liquid)*Item No.: MED23786*

Purity Specification:	Plasma Fraction
Source:	Pooled Human Donors
Total Protein:	15.0 - 40.0 g/dL (Biuret)
Appearance:	Yellow to orange liquid
Total Cholesterol:	Typically undetectable
Triglycerides:	Typically undetectable
Storage:	-20°C , 2 years

Infectious Disease Testing:

Negative or non-reactive at the donor level for HIV-1, HIV-2 and HCV (antibodies), HIV-1, HBV and HCV (NAT), HBsAg and syphilis according to FDA approved methods.

Conclusion: This batch of products conform the quality inspection standard of the manufacturer