

Product Name: D-Dimer antigen(Liquid)*Item No.: MED23463*

Purity Specification:	95% or 50% (SDS-PAGE)
Source:	Human Plasma
D-Dimer Content:	1 - 5 mg/mL (ELISA)
Appearance:	Clear and colorless solution
Molecular Weight:	~200,000
Storage:	-20°C , 2 years

Formulation:

0.2 micron filtered solution in 0.05 M sodium phosphate, 0.15 M sodium chloride, pH 7.5 with 0.09% sodium azide

Infectious Disease Testing:

Negative or non-reactive at the donor level for anti-HIV 1 and 2, anti-HCV, HBsAg, HIV-1 (NAT), HBV (NAT), HCV (NAT), and syphilis by FDA approved methods.

Conclusion: This batch of products conform the quality inspection standard of the manufacturer