

Product Name: NT-proBNP (27-31) Antibody*Item No.: MEM24399*

Purity Specification:	90%
Source:	Mouse
Format:	Purified Monoclonal Antibody
Immunogen:	Native Human ProBNP antigen
Applications:	Lateral Flow
Supplied In:	10mM Phosphate Buffered Saline, pH 7.4
Storage:	Short Term (≤ 2 weeks): 2-8°C Long Term: -20°C Avoid repeated freezing and thawing

Conclusion: This batch of products conform the quality inspection standard of the manufacturer