

Sample Name: Fusaric Acid

CAS NO: 536-69-6

Acq. Operator: Armay Chen

Acq. Instrument: Instrument 12

Date Acquired: 1/19/2022 5:51:04 PM CST

Vial: 2:B,2

Acq. Method Set: 1_POS_2MIN

Date Processed: 1/20/2022 3:52:29 PM CST

Injection Volume: 0.50 ul

Processing Method: PDA_01, MS_scan




Run Time: 2.8 Minutes

Channel Name: PDA Ch1 214nm@4.8nm, MS TIC

Column Name: Poroshell 120 EC-C18 A-RP-929

Injection #: 1

Area Percent Report



Peak Results

Channel Name: PDA Ch2 254nm@4.8nm

	RT	Width (sec)	Area	Height	% Area	Channel Name
1	0.981	1.800	2199	2505	0.42	PDA Ch2 254nm@4.8nm
2	1.022	6.700	522803	404291	99.23	PDA Ch2 254nm@4.8nm
3	1.675	2.300	1845	1600	0.35	PDA Ch2 254nm@4.8nm

Peak Results

Channel Name: PDA Ch1 214nm@4.8nm


	RT	Width (sec)	Area	Height	% Area	Channel Name
1	0.981	1.850	4643	4765	0.55	PDA Ch1 214nm@4.8nm
2	1.022	6.100	826017	643353	98.63	PDA Ch1 214nm@4.8nm
3	1.603	2.000	4591	4339	0.55	PDA Ch1 214nm@4.8nm
4	1.673	2.500	2223	1751	0.27	PDA Ch1 214nm@4.8nm

Peak Results

Channel Name: MS TIC

	RT	Width (sec)	Area	Height	% Area	Channel Name
1	1.097	15.958	583193613	105367871	100.00	MS TIC

Match Plot



*** End of Report ***