


=====
Acq. Operator : Armay Chen
Acq. Instrument : Instrument 1100
Injection Date : 11/29/2021 10:39:43 AM
Acq. Method : C:\CHEM32\MED12081\N1120A.M
Last changed : 11/29/2021 12:44:56 AM by Armay Chen
(modified after loading)
Analysis Method : C:\CHEM32\1\METHODS\GENERIC.M

Sample Info : Colivelin (trifluoroacetate salt); Cat #MED12081; QC
Column:Unitary C18, 4.6*150mm, 5um 100A
Buffer A:0.1%TFA in H2O; Buffer B:0.1%TFA in H2O
Gradient:31-41% B in20min , Flow rate:1ml/min



=====
Area Percent Report
=====

Signal 1: DAD1 A, Sig=220, 4 Ref=360, 100

Peak #	RetTime [min]	Type	Width [min]	Area [mAU*s]	Height [mAU]	Area %
1	6.989	VV	0.3520	24.3241	0.9096	1.0172
2	8.654	VV	0.2256	29.7513	1.6250	1.2442
3	9.460	VV	0.3817	2287.1230	90.8362	95.6454
4	10.439	VB	0.3835	50.0550	1.6806	2.0933

Totals : 2391.2534

=====
*** End of Report ***
=====