

Sample Name: TMS

Acq. Operator: Armay Chen

Acq. Instrument: HPLC-1260

Vial: 1:D,6


Run Time: 10.0 Minutes

Channel Name: 254nm

Inj volume: 0.20 ul

Proc.Chnl.Descr: PDA Spectrum PDA 254.0nm

### Area Percent Report



\_\_\_ Processed Channel Descr. PDA Spectrum PDA 254.0 nm (PDA Spectrum (210-400)nm)

| Rank | RT    | Area   | Height | Area % |
|------|-------|--------|--------|--------|
| 1    | 3.705 | 9998   | 3267   | 1.58   |
| 2    | 6.726 | 619584 | 237058 | 98.10  |
| 3    | 7.009 | 1989   | 1030   | 0.31   |

\*\*\* End of Report \*\*\*