

Data File C:\CHEM32\MED18397\M2G8397CA.D

Sample Name: Aprotinin

Acq. Operator: Armay Chen

Acq. Instrument: Instrument 11

Column: Venusil ASB C18, 4.6*250mm, 5 μ m, 300A

Solvent A: 0.1%Trifluoroacetic in 100% Acetonitrile

Solvent B: 0.1%Trifluoroacetic in 100% Water

Flow rate: 1.0ml/min


Wavelength: 220nm

Volume: 10 μ l

Gradient:

Time	A	B
0.0min	5%	95%
30.0min	95%	5%
35.0min	Stop	Stop

Area Percent Report



Rank	RetTime(min)	Area(mAU*S)	Area %
1	11.651	654983	99.304
2	14.499	5412	0.696

*** End of Report ***